



Beauty Teck

International (Pvt) Ltd.



**Catalogue and Key performance claims of
Non-Sterile, Single Use, Retractors**

General Instruments Retractors



Key performance claims of Non-Sterile, Single Use, Retractors

1. Our Retractors are provided good exposure without the need of multiple surgical assistants.
2. Device is Biocompatible
3. Device is free of sharp edges, pits and deep scratches
4. Material is Stainless Steel AISI 410
5. Device bears acceptable hardness (40-45 HRC)
6. Corrosion Free Device
7. Shelf life is 3 years



Weitlaner
Sharp
10.5 cm/4¼"
17-05781-10



Weitlaner
Sharp
14.0 cm/5½"
17-05781-14



Weitlaner
Blunt
10.5 cm/4¼"
17-05782-10



Weitlaner
Blunt
14.0 cm/5½"
17-05782-14



Weitlaner
Sharp
14.0 cm/5½"
17-05785-14




Weitlaner
Sharp
16.5 cm/6½"
17-05785-16

General Instruments Retractors



Key performance claims of Non-Sterile, Single Use, Retractors

1. Our Retractors are provided good exposure without the need of multiple surgical assistants.
2. Device is Biocompatible
3. Device is free of sharp edges, pits and deep scratches
4. Material is Stainless Steel AISI 410
5. Device bears acceptable hardness (40-45 HRC)
6. Corrosion Free Device
7. Shelf life is 3 years



Weitlaner
Blunt
14.0 cm/5½"
17-05786-14



Weitlaner
Blunt
16.5 cm/6½"
17-05786-16



Weitlaner
Blunt Retractor
17-05793-15
16 cm

Beauty Teck

International (Pvt) Ltd.

Head Marala Road, Machi Khokhar
Sialkot, 51310-Pakistan.

Tel: 0092 52 4295500, 4564343

Fax: 0092 52 4298551

Email: info@beauty-teck.com

Website: www.beauty-teck.com