



# Beauty Teck

International (Pvt) Ltd.

Brochure and Key performance claims of Sterile & Non-Sterile  
**Single Use Scalpel Handles**







# SCALPEL HANDELS

## Key performance claims of Sterile & Non-Sterile, Single Use, Scalpel Handles

- Claim 1 Our Scalpel Handles fitment is accurate with relevant blades
- Claim 2 Device is Biocompatible
- Claim 3 Device is free of sharp edges, pits and deep scratches
- Claim 4 Material is Stainless Steel AISI 410
- Claim 5 Device bears acceptable hardness (40-45 HRC)
- Claim 6 Corrosion Free Device
- Claim 7 Shelf life is 5 years for Sterile and 2 years for non-sterile items
- Claim 8 ETO Sterilized Product



**No 3**  
**13-01348-03**




# SCALPEL HANDELS

## Key performance claims of Sterile & Non-Sterile, Single Use, Scalpel Handles


- Claim 1 Our Scalpel Handles fitment is accurate with relevant blades
- Claim 2 Device is Biocompatible
- Claim 3 Device is free of sharp edges, pits and deep scratches
- Claim 4 Material is Stainless Steel AISI 410
- Claim 5 Device bears acceptable hardness (40-45 HRC)
- Claim 6 Corrosion Free Device
- Claim 7 Shelf life is 5 years for Sterile and 2 years for non-sterile items
- Claim 8 ETO Sterilized Product




**FRONT VIEW**



**SIDE VIEW**



**Enlarged View**




**No 3**  
**13-01352-03**




# SCALPEL HANDELS

## Key performance claims of Sterile & Non-Sterile, Single Use, Scalpel Handles


- Claim 1 Our Scalpel Handles fitment is accurate with relevant blades
- Claim 2 Device is Biocompatible
- Claim 3 Device is free of sharp edges, pits and deep scratches
- Claim 4 Material is Stainless Steel AISI 410
- Claim 5 Device bears acceptable hardness (40-45 HRC)
- Claim 6 Corrosion Free Device
- Claim 7 Shelf life is 5 years for Sterile and 2 years for non-sterile items
- Claim 8 ETO Sterilized Product




**ENLARGED VIEW**



**FRONT VIEW**



**SIDE VIEW**



**No 4**  
**13-01356-04**



# Beauty Teck

International (Pvt) Ltd.

☎ +92-52-4564343, 4363030

🌐 [www.beauty-teck.com](http://www.beauty-teck.com)

☎ +92-321-6101111

✉ [info@beauty-teck.com](mailto:info@beauty-teck.com)

📍 Head Marala Road, Machi Khokhar Sialkot-51610-Pakistan